

MRI Accreditation Program Visual Checklist

MRI Facility Name:

MRI Scanner Identifier:

						Date:
	Table position and other displays					
Patient Transport and Gantry	Alignment lights					
	Horizontal smoothness of motion and stability					
	Vertical motion smoothness and stability					
Filming Viewing	Laser camera					
	Light boxes					
RF Integrity and Control Room	RF door contacts					
	RF window-screen integrity					
	Operator console switches/lights/meters					
	Patient monitor (if present)					
	Patient intercom					
	Room temperature/room humidity					
	Emergency cart					
	Safety warning signage					
Facility Safety	Door indicator switch (if installed)					
	Cryogen level indicator					
Pass = ☒ Fail = F Does Not Apply = NA						
Technologist Initials:						

Reviewed by:

Qualified Medical Physicist/MRI Scientist

Date of Review:

WEEKLY MRI EQUIPMENT QUALITY CONTROL FOR LARGE PHANTOM

MRI Facility Name:

MRI Scanner Identifier:

[illegible]

Reviewed by:

Qualified Medical Physicist/MRI Scientist

Date of Review: